

PURAXEL

Laser

• Ultrasonic

• Galvanic

• Cooling

• LED



2940nm Er:YAG
Fractional Laser
Aesthetic beauty device



Laser beauty device
makes your own way for skincare
www.puraxel.co.kr



All · in · One

“ No needle!
No pain!
No 2nd infection! ”



* This device is a skin care device, not a medical device.

LASER BEAUTY DEVICE

MAKES YOUR OWN WAY FOR SKINCARE

- Laser Type : Er:YAG
- Laser Class : Class 1
- Laser Level : Steps 1~3
- Weight : Main Body 520g / Cradle 410g
- Dimensions : Body 82(W)x201(D)x225.3(H)mm / Cradle 250x169x135mm
- Battery : 11.1V 1,300mAh Lithium-Polymer Battery
- Charging Adapter : 14VDC / 1.71A

Laser Mode

level: 1~3

Creates micro-pores in the skin
Create multi micro holes in the skin and dramatically deliver nutrients and various cosmetics.

Ultrasonic Mode

treatment 5 minutes

Functional nutrient absorption.
Improve absorption rate of beauty solutions by ultrasonic waves.

Galvanic Mode Option

treatment 4 minutes

Deep cleansing and nutrient absorption
Positive(+) nature of the rectified current, discharges waste, while the negative(-) nature enhances the delivery rate of functional nutrients.

Cooling Mode Option

treatment 3 minutes

Soothing skin and reducing pores
Apply cooling mode to minimize and settle pores in the skin after using laser, ultrasonic, or LED mode.

LED Mode Option



treatment 3 minutes

Moisturizing, elasticity and soothing effect
Through the Red, Blue and Pink modes, calm troubled skin (acne), moisturize, and improve elasticity

CLINICAL LAB RESULTS

+ Moisture / Density / Tone Brightness

* The efficacy of laser beauty device with a cosmetic product




1 Three times the increase in Skin Moisture Content.

2 30% increase in Skin Density.

3 2 times improvement in Skin Tone Brightness.

+ Melasma / Wrinkles / Hair loss



Melasma improvement

Wrinkles improvement

Hair loss improvement



LA MEDITECH #1002, ACE High-End 6th, Beotkkot-ro 234, Geumcheon-gu, Seoul, South Korea 08513 · Tel : +82-2-852-7982 · Fax : +82-2-852-7983

www.lameditech.com · E-mail : puraxel@lameditech.com